

## ????????????????? ?????????? KERN ORA-P



????????????????? ?????????? KERN ORA-P Analogue refractometer ORA-P Area of application: Urine

?????????????: ??? ?????????????? ??????

????

???? ?? ??????? 80,00 €

???? ??????? 99,20 €

???? ??????? ?????? 80,00 €

????????

????????? ?????? 19,20 €

[????????? ?? ?? ??????](#)

??????????

### ????????????????? ?????????? KERN ORA-P

Analogue refractometer ORA-P Area of application: Urine

- The models in the KERN ORA-P range are universal, maintenance-free analogue manual refractometers
- The practical and robust design means that it can be used simply, efficiently and on a permanent basis, day in day out
- The manual effort to convert the value is avoided by the use of several selectable scales and thereby eliminates user error
- These scales have been specially developed, calculated precisely and checked. Another key feature is the very thin, clear lines
- The optical system and the prism cover are made from special materials, which enable you to measure with low tolerances
- All models are fitted with an eyepiece which can be adjusted easily and smoothly to different visual capacities
- The models which are marked "ATC" have automatic temperature compensation, which makes it possible to obtain precise measurements at different environmental temperatures ( °C - °C)
- Included in the scope of supply:
  - Storage box
  - Calibration solution
  - if necessary, calibration block (only for KERN ORA 2PB)
  - Pipette
  - Screwdriver
  - Cleaning cloth
- Further accessories are available as options

Area of application: Urine

The following models are particularly suitable for measuring specific urine weight (density), the serum content (serum protein in urine) and the refraction index.

Main areas of application:

- Hospitals
- Doctor's surgeries
- Medical training institutions
- Care homes and nursing homes
- Sports medicine (drug testing)

?? ?????????? ?????????????? **ORA-P** ?????? ?????????? ?? ?????????? ?????? ??? ?? ???????????? ??????????????:

[?????????????.....](#)